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iMethod™ test kit

- method documentation
- MS and LC method files
- processing method files and report templates
- sample preparation
- instrument parameters (to accelerate the adoption of this method for routine testing)

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LC

- System controller CBM-20A
- 2 Isocratic pumps LC-20AD (with semi-micro 50 µL mixer)
- Autosampler SIL-20AC
- Column oven CTO-20AC

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

- 1. Measure a volume of urine sample into a microcentrifuge tube.
- 2. Add x μ L of combined internal standard spiking solution to each tube and mix using a vortex mixer for approximately 5 to 10 seconds.
- 3. Add x μ L of methanol to each double blank tube and mix using a vortex mixer for approximately 5 to 10 seconds.
- 4. Add x μ L of water to each tube and mix using a vortex mixer for approximately 5 to 10 seconds.
- 5. Condition the SPE phase (cartridge or 96 well plate) with x mL of SPE Wash solution.
- 6. Apply vacuum or positive pressure.
- 7. Condition the SPE phase (cartridge or 96 well plate) with x mL of water.
- 8. Apply vacuum or positive pressure.
- 9. Add sample and apply vacuum or positive pressure.
- 10. Wash the SPE phase (cartridge or 96 well plate) with x mL of water.
- 11. Add sample and apply vacuum or positive pressure.
- 12. Wash the SPE phase (cartridge or 96 well plate) with 2 x 1 mL of SPE Wash solution.
- 13. Apply vacuum or positive pressure for each 1 mL addition.
- 14. Replace tubes or 96 well plate in preparation for elution.
- 15. Elute from the SPE phase (cartridge or 96 well plate) with 3 x 0.35 mL of SPE Elution

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


Data File C:\HPCHEM\1\DATA\260412\MET15_04.D Sample Name: CS 1810

Injection Date : 4/26/2012 6:20:20 PM Seq. Line : 3
Sample Name : CS 1810 Vial : 41
Acq. Operator : yc Inj : 1
Inj Volume : 40 µl
Sequence File : C:\HPCHEM\1\SEQUENCE\META15YC.S
Method : C:\HPCHEM\1\METHODS\META.M
Last changed : 10/14/2011 1:54:58 PM by Yc

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