

Analysis of Dioxins and Dioxin-like PCBs in Feed and Food

Thierry Faye, Agilent Technologies
Jef Focant, University of Liege, Belgium

GC/MS/MS Analyzer for Analysis of Dioxins and Dioxin-like PCBs in Feed and Food

Thierry Faye
Market Development Manager
Agilent Technologies

Agilent Technologies

NEW Introducing the Agilent GC/MS/MS Dioxin Analyzer

→ Cutting edge technology enables fast & reliable Dioxin analysis at low levels:

- New standard for GC/MS/MS sensitivity with 7010 Triple Quadrupole EI source = An instrument with sensitivity of GC Sector MS.
- MultiMode Inlet (MMI) for effective cold splitless injections and more.

→ Developed from successful collaborations with leading Dioxin Labs in Europe:

- The Agilent platform already validated according to new regulations in Europe for both food and feed (EC 589/2014, 709/2014).
- Refer to L'Homme/Focant's publication: *Journal of Chromatography A*, 1376 (2015) 149-158.
- Custom reporting with complete calculations have been developed and automated in MassHunter.

→ Ready for analysis:


- Pre-configured and pre-tested at our factory so installation in your lab is fast and efficient.
- The RTL advantage:** guarantees the exact matching of our reference method on a new instrument.
- A service engineer runs a complete check out standard so validation can begin.
- Method never needs altering even when column maintenance is performed.

Agilent Technologies | 2-15-17

Validation of GC/MS/MS confirmatory method for the European official control of levels of dioxins, furans, and dioxin-like PCBs in foodstuffs

L'Homme B., Focant J.-F.

Organic and Biological Anal. Chem.
University of Liege, Belgium



95% Exposure by Food Consumption

➔ **Food and Feed control...**



Dioxin EU Regulation for food-feed 'started' with the Belgian Dioxin Crisis in 1999...

EU Commission Documents

- Council Regulation 2001/102/EC
- Council Regulation (EC) No 2375/2001
- Commission Directive 2002/69/EC
- Commission Directive 2002/70/EC
- Commission Directive 2006/13/EC
- Commission Recommendation 2006/181/EC
- Commission Regulation (EC) No 1881/2006
- Commission Regulation (EC) No 1883/2006
- Commission Regulation (EC) No 152/2009
- Commission Regulation (EU) No 1259/2011
- Commission Recommendation 2011/516/EU
- Commission Regulation (EC) No 252/2012
- Commission Regulation (EC) No 277/2012
- Commission Regulation (EC) No 278/2012
- **Commission Regulation (EC) No 589/2014**
- Commission Regulation (EC) No 709/2014



<http://eur-lex.europa.eu>

EU Commission Strategy (food-feed)

- ✓ Continuous monitoring
- ✓ Maximum-Action-(Target) level strategy
- ✓ Screening-Confirmatory approach
- ✓ RASFF (high capacities)

➔ Based on state-of-the-art methods

Analytical Methods ?

Evolutionary Guidelines

HARMONISED QUALITY CRITERIA FOR CHEMICAL AND BIOASSAYS ANALYSES OF PCDDs/PCDFs IN FEED AND FOOD PART 1: GENERAL CONSIDERATIONS, GC/MS METHODS

Rainer Malsch¹, Bert Baumann¹, Peter A. Behnisch¹, Richard Canady², Daniel Fraisse¹, Peter Fürst¹, Douglas Hayward², Ronald Hoogenboom¹, Ronald Hoogerbrugge², Eijun Imai³, Otho Pappas⁴, Wim Traug¹, Thomas Wilmüller¹

OHCE 50 (2001) 53

➔ **PBMS**

↓

HARMONISED QUALITY CRITERIA FOR CHEMICAL AND BIOASSAYS ANALYSES OF PCDDs/PCDFs IN FEED AND FOOD PART 2: GENERAL CONSIDERATIONS, BIOASSAY METHODS

Peter A. Behnisch¹, Randy Allen², Jack Anderson³, Abraham Brouwer¹, David J. Brown⁴, T. Colin Campbell⁵, Leo Goeyens¹, Robert O. Harrison¹, Ron Hoogenboom¹, Ilse Van Overmeire¹, Wim Traug¹, and Rainer Malsch¹

OHCE 50 (2001) 59

- GC-IDHRMS vs CALUX
- ISO 17025
- Validation @ LOQs (@1/5th level of interest)
- #Upper-Lowerbound < 20%
- @1pg TEQ/g fat level
- Z-scores @ PTs
- Recovery rates...

www.dioxin20xx.org

Recent EN 16215

NEN-EN 16215:2012

EUROPEAN STANDARD EN 16215

NORME EUROPÉENNE

EUROPÄISCHE NORM April 2012

ICS 65.120


English Version

Animal feeding stuffs - Determination of dioxins and dioxin-like PCBs by GC/HRMS and of indicator PCBs by GC/HRMS

Aliments des animaux - Dosage des dioxines, des PCB de type dioxine et des PCB indicateurs par GC/HRMS Futtermittel - Bestimmung von Dioxinen und dioxinähnlichen PCBs mittels GC/HRMS und von Indikator-PCBs mittels GC/HRMS

This European Standard was approved by CEN on 9 March 2012.

PBMS philosophy



Confirmatory Tool

✓ HRMS sector or MS/MS instruments

ANALYTICAL CRITERIA FOR USE OF MS/MS FOR DETERMINATION OF DIOXINS AND DIOXIN-LIKE PCBs IN FEED AND FOOD

EURL Working group recommendations

Kotz A¹, Malisch R^{1†}, Focant J², Eppe G², Cederberg TL³, Rantakokko P⁴, Finnsson O⁵, G. J. Van Leeuwen J⁶, L. Lovász C⁷, Scorticchi G⁸, Dieltz G⁹, di Domenico A¹⁰, Inesido A¹¹, K. M. M. van Leeuwen J¹², M. A. M. van Leeuwen J¹³

¹European Union Reference Laboratory (EU-RL) for Dioxins and PCBs, National Institute for Research in Food, Freiburg, Germany; ²CART, University of Liege, Belgium; ³National Institute for Environmental Health and Safety, Copenhagen, Denmark; ⁴National Institute for Environmental Health and Safety, Helsinki, Finland; ⁵Chemisches und Veterinäruntersuchungsamt Münster-Lippe (CVUA-MEL), Münster, Germany; ⁶NCSR Department of Analytical Chemistry, Institute for Agricultural and Fisheries Research, Teagasc, Dublin, Ireland; ⁷Central Institute for Food Safety and Food Quality, Hungary; ⁸Istituto Zooprofilattico Sperimentale dell’Abruzzo e Molise “G. Caporale”, Teramo, Italy; ⁹Istituto Superiore di Sanità (ISS), Rome, Italy; ¹⁰Central Laboratory for Food Safety and Food Quality, Wageningen, The Netherlands; ¹¹The Food and Environment Research Agency (FERA), York, United Kingdom

Organohalogen compound 74 (2012), 156-159

EU Regulation 252/2012 amended by **589/2014** (2nd of June 2014)
 EU Regulation 152/2009 amended by **709/2014** (20th of June 2014)

New EU Regulation

L 164/18 Official Journal of the European Union 3.6.2014

COMMISSION REGULATION (EU) No 589/2014
of 2 June 2014

laying down methods of sampling and analysis for the control of levels of dioxins, dioxin-like PCBs and non-dioxin-like PCBs in certain foodstuffs and repealing Regulation (EU) No 252/2012
(Text with EEA relevance)

27.6.2014 Official Journal of the European Union L 188/1



COMMISSION REGULATION (EU) No 709/2014
of 20 June 2014

amending Regulation (EC) No 152/2009 as regards the determination of the levels of dioxins and polychlorinated biphenyls
(Text with EEA relevance)


(9) In addition to the gas chromatography/high resolution mass spectrometry (GC-HRMS), technical progress and developments have shown that also gas chromatography/tandem mass spectrometry (GC-MS/MS) can be used as a confirmatory method for checking compliance with the maximum level (ML). Regulation (EU) No 252/2012 should therefore be replaced by a new Regulation providing for the use of gas chromatography/tandem mass spectrometry (GC-MS/MS) as an appropriate confirmatory method for checking compliance with the maximum level.

| Criteria | PCDD/Fs and DLPCBs GC-MS/MS (589/2014) | NDLPCBs GC-MS/MS (709/2014) |
|-------------------------|--|---|
| Detectable quantity | PCDD/F upper femogram (10 ⁻¹⁵ g) NDLPCB low picogram (10 ⁻¹² g) MDLPCB nanogram (10 ⁻⁹ g) | NDLPCB nanogram (10 ⁻⁹ g) |
| Selectivity | Chromatographic separation of 2,3,4,7,8-HxCDF and 1,2,3,6,8-HxCDF ≥5% valley peak to peak | Relative RT ±0.25% 15-s analysis |
| MRM transitions | Monitoring 2 specific precursors with both specific product ion transition for all labeled and unlabeled analysis Relative ion intensity max ±15% Resolution MS quadrupole = unit | Monitoring at least 1 precursor ion and 2 product ions Tolerance ratio ±20% if rel. intens. ≥0.6 Tolerance ratio ±25% if rel. intens. 0.20-0.6 Resolution MS quadrupole = unit |
| Blank | Used for LOD calculation | Used for LOD calculation |
| LOD | LOD calculated from lowest cal. point lowest concentration point on cal. multiple applicable and consistent relation to the average RRF Average RRF calculated for all points Deviation to average RRF < 0% | ditto |
| LOQ | LOQ calculated from average blank level LOQ = 1.0 of maximum ML Difference ub and lb loads < 0.5 ML Demonstrate performance at 0.5 ML ML 2 ML | ditto |
| Accuracy | Recovery accuracy ±20% | Demonstrate performance at 0.5 ML ML 2 ML |
| Reproducibility | Within-lab reproducibility (RSD) < 5% | Recovery for sumind-PCB @ ML ±20% Within-lab reproducibility (RSD) < 2% |
| Control | QC charts for blanks QC charts control sample | QC charts for blanks QC charts control sample |
| Recovery | Individual internal std in range 40-120% Out of range OK if contrib. to 100 ± 0% | Individual internal std in range 50-120% Out of range OK if contrib. to sumind-PCB 0% |
| Measurement uncertainty | Coverage factor = 2 (C165%) If separate determination of components, value = sum of separate uncertainty for sum of PCDD/F and DLPCBs | Expanded measurement uncertainty Coverage factor = 2 (C165%) |

From Regulation
↓
Full validation




Chromatographic Profile



Chromatographic Separation


Selectivity



✓ 25% valley separation HxCDF

Tandem in-space MS

7000 MS/MS System
Optimized for gas chromatography



| PCDD/Fs | PCB | PCDD/Fs | PCB | PCDD/Fs | PCB | PCDD/Fs | PCB |
|--------------|-------|-------------|-------|-------------|-------|-------------|-------|
| 12378-TCDD | 25.72 | 22478-PeCDF | 24.08 | 12378-TCDF | 26.20 | 12378-PeCDF | 23.24 |
| 12378-PeCDD | 24.33 | 12378-PeCDF | 27.02 | 12378-PeCDF | 27.51 | 22478-PeCDF | 27.70 |
| 12378-TCDF | 27.93 | 12378-PeCDF | 28.64 | 12378-PeCDF | 31.95 | OCDF | 31.62 |
| 12378-PeCDF | 26.49 | OCDF | 35.57 | OCDF | 35.97 | | |
| 123878-TCDD | 32.26 | | | | | | |
| 123878-PeCDF | 32.26 | | | | | | |
| 12378-TCDF | 26.20 | | | | | | |
| 12378-PeCDF | 23.24 | | | | | | |
| 22478-PeCDF | 24.08 | | | | | | |
| 12378-PeCDF | 27.02 | | | | | | |
| 123878-PeCDF | 27.51 | | | | | | |
| 22478-PeCDF | 27.70 | | | | | | |
| 12378-PeCDF | 28.64 | | | | | | |
| 123878-PeCDF | 31.95 | | | | | | |
| OCDF | 31.62 | | | | | | |
| OCDF | 35.76 | | | | | | |

MRM transition Ratio

Quant/Qual transition ratio

EU Reg 709/2014 says:

✓ PCDD/Fs, DL-PCBs: 2 specific precursors with each specific product ions

• 2378 TCDD: 319.9 > 256.9 → Average ratio = 0.964
321.9 > 258.9

✓ NDL-PCBs: at least 1 precursor and 2 product ions

• PCB 189: 393.8 > 323.8 → Average ratio = 0.627
395.8 > 325.8

Tolerance

± 15%
R(quad)=unit

± 20%
R(quad)=unit

In MassHunter



Proper Estimation of LOQs

LOD/LOQ

ANALYTICAL CRITERIA FOR USE OF MS/MS FOR DETERMINATION OF DIOXINS AND DIOXIN-LIKE PCBs IN FEED AND FOOD

Kotz A¹, Maljoch R², Focant J², Eppe G², Cedergren TL³, Rannakko P⁴, First P⁵, Bernbaum T⁴, Leonhardt L⁴, Lovisa C⁶, Scottichini G⁶, Dieltz G⁶, Di Domenico A⁶, Ingelido AM⁷, Traag W⁸, Smith F⁴, Fernandez A¹¹

Approach 1: The LOQ can be calculated from the signal-to-noise ratio as already defined in the current regulations.

Approach 2: As an alternative approach, if the signal-to-noise ratio does not provide reliable results due to a very low, or no discernable noise level, the limit of quantification is based on the calibration curve. The limit of quantification is then defined as the lowest concentration point on a calibration curve that gives an acceptable (< 30 %) and consistent (measured at least at the start and at the end of an analytical series of samples) deviation to the average relative response factor calculated for all points on the calibration curve in each series of samples.

For conversion of the limit of quantification from the calibration curve to the sample, the recovery of the internal standards of the respective congener and the sample intake has to be taken into account.

Organohalogen compound 74 (2012), 156-159

- ✓ iLOQ = 10*stdev (8 replicate injections – cali point)
- ✓ LOQ = blank mean + 6*stdev (12 injections – blank)

Validation for vegetable oil (feed)

❖ iLOD/iLOQ

'Acceptable and consistent deviation to the average RRF'

| Quantifier | Name | TS | RT | Transition | Scan | Type |
|-------------|-------|--------|----------|------------|--------|-----------|
| Calibration | | | | | | |
| | Level | Conc. | Response | Enable | RF | Final RSD |
| 1 | CDC01 | 0.0160 | 424 | ✓ | 0.9341 | 2.534717 |
| 2 | CDC01 | 0.0160 | 450 | ✓ | 0.9307 | 2.524717 |
| 3 | CDC01 | 0.0160 | 460 | ✓ | 0.9257 | 2.534717 |
| 4 | CDC02 | 0.0400 | 398 | ✓ | 0.9298 | 2.534717 |
| 5 | CDC02 | 0.0400 | 844 | ✓ | 1.0297 | 7.235217 |
| 6 | CDC03 | 0.0800 | 2710 | ✓ | 0.9728 | 3.487891 |
| | CDC03 | 0.0800 | 2873 | ✓ | 0.9276 | 3.487891 |
| | CDC04 | 0.4000 | 13004 | ✓ | 0.9468 | 2.407250 |
| | CDC04 | 0.4000 | 1304 | ✓ | 0.9772 | 2.407250 |
| | CDC05 | 0.8000 | 26258 | ✓ | 0.9742 | 1.207790 |
| | CDC05 | 0.8000 | 26255 | ✓ | 0.9670 | 1.207790 |
| | CDC05 | 0.8000 | 26843 | ✓ | 0.9801 | 1.207790 |
| | CDC06 | 2.0000 | 89292 | ✓ | 0.9997 | 3.687836 |
| | CDC06 | 2.0000 | 89441 | ✓ | 1.0411 | 3.682236 |
| | CDC06 | 2.0000 | 89721 | ✓ | 1.0295 | 3.682236 |

Av RRF₁₋₆ = 0.9622
 $R^2(\text{linear fit}) = 0.996$
 > 0.9900
 Av RRF_i = 0.9435
 Deviation = -1.9%
 $< 30\%$

Validation for vegetable oil (feed)

❖ iLOD/iLOQ

Calibration curve (lowest level)

8 replicates

10*stdev

❖ LOQ matrix

Procedure blanks

10 independent injections

Average + 6*stdev.

Validation for vegetable oil (feed)

❖ (i)LOQ

| Compound | Lowest call point | | | | Blanks | | |
|--------------|-------------------|-----------|-----------|-----------|----------|----------------|---------------|
| | Avg Conc pg/g | Std. Dev. | iLOD pg/g | iLOQ pg/g | In blank | iLOQ MS/MS PpB | iLOQ HRMS PpB |
| PCB 77 | 0.2917 | 0.0037 | 0.011 | 0.037 | yes | 49.66 | 64.59 |
| PCB 81 | 0.2935 | 0.003 | 0.009 | 0.030 | yes | 3.67 | 9.53 |
| 12378-TCDF | 0.0152 | 0.0018 | 0.005 | 0.018 | no | 0.02 | 0.06 |
| 12378-TCDF | 0.0162 | 0.001 | 0.003 | 0.010 | yes | 0.10 | 0.12 |
| PCB 126 | 0.2906 | 0.0077 | 0.023 | 0.077 | yes | 1.37 | 1.21 |
| 123478-PCDF | 0.0154 | 0.0021 | 0.006 | 0.021 | yes | 0.08 | 0.08 |
| 12378-PCDF | 0.0169 | 0.0029 | 0.009 | 0.029 | no | 0.01 | 0.06 |
| PCB 169 | 0.3062 | 0.0071 | 0.021 | 0.071 | no | 0.02 | 0.27 |
| 12378-PCDF | 0.0142 | 0.0022 | 0.007 | 0.022 | yes | 0.56 | 0.09 |
| 123478-PCDF | 0.0164 | 0.0016 | 0.005 | 0.016 | yes | 0.07 | 0.10 |
| 123678-PCDF | 0.0169 | 0.0009 | 0.003 | 0.009 | yes | 0.06 | 0.06 |
| 1234678-PCDF | 0.0166 | 0.0007 | 0.002 | 0.007 | yes | 0.08 | 0.09 |
| 123789-PCDF | 0.0156 | 0.002 | 0.006 | 0.020 | yes | 0.18 | 0.10 |
| 123789-PCDF | 0.0174 | 0.0061 | 0.018 | 0.062 | yes | 0.06 | 0.06 |
| 123478-PCDF | 0.0117 | 0.0022 | 0.007 | 0.022 | yes | 0.03 | 0.06 |
| 123678-PCDF | 0.0426 | 0.0032 | 0.010 | 0.032 | yes | 0.14 | 0.08 |
| PCDF | 3.5201 | 0.0465 | 0.140 | 0.465 | yes | 1.99 | 2.02 |
| PCDF | 0.0158 | 0.0027 | 0.008 | 0.027 | yes | 0.60 | 0.77 |

Accuracy & Reproducibility

Spiked materials at 0.5 ML, ML, 2 ML (6 series, 3 days)


| PCDD/Fs | Spike level | Average | Stdev | RSD | Target | Bias |
|---------|-------------|---------------|---------------|-----|---------------|-------|
| | | ng WHO-TEQ/kg | ng WHO-TEQ/kg | % | ng WHO-TEQ/kg | % |
| ML/2 | ML | 0,409 | 0,029 | 7,1 | 0,40 | 2,36 |
| | ML | 0,778 | 0,045 | 5,7 | 0,79 | -1,54 |
| | 2ML | 1,600 | 0,035 | 2,2 | 1,58 | 1,30 |
| NO-PCBs | Spike level | Average | Stdev | RSD | Target | Bias |
| | | ng WHO-TEQ/kg | ng WHO-TEQ/kg | % | ng WHO-TEQ/kg | % |
| ML/2 | ML | 0,307 | 0,028 | 9,0 | 0,33 | -7,00 |
| | ML | 0,595 | 0,020 | 3,4 | 0,65 | -8,53 |
| | 2ML | 1,256 | 0,021 | 1,6 | 1,30 | -3,42 |

- ✓ Bias < 20%
- ✓ Within lab reproducibility < 15%

Accuracy & Reproducibility


Proficiency test (PT)

| Material | Reported pg/g TEQ | Target Value pg/g TEQ | Accuracy |
|-------------------|-------------------|-----------------------|----------|
| Material 1 | | | |
| PCDD/Fs | 1.10±0.20 | 1.01 | 8.8% |
| DL-PCBs | 0.80±0.19 | 0.89 | -10.3% |
| Total TEQ | 1.90±0.36 | 1.90 | -0.1% |
| Material 2 | | | |
| PCDD/Fs | 0.55±0.11 | 0.48 | 15.6% |
| DL-PCBs | 0.82±0.21 | 0.85 | -3.0% |
| Total TEQ | 1.38±0.26 | 1.33 | 3.7% |



Within Lab Reproducibility

QC pork fat



MassHunter Report Generator

| D:\MassHunter\GCMS1\Data\Dev\pCDDF-NO\PCDD\Scan Results\LabData\PCDD-F-NO1.mh.ms | | | | | | | | | |
|--|------------------|---------|------------------|-------------|------------------|----------|------------------|---------------|--|
| Batch Data Path | 2014-06-12 10:33 | Analyst | admin | Report Time | 2014-06-12 13:30 | Reporter | admin | Lab Data Path | D:\MassHunter\GCMS1\Data\Dev\pCDDF-NO\PCDD\Scan Results\LabData\PCDD-F-NO1.mh.ms |
| Law Calif Update | 2014-06-12 10:33 | Batch | 2014-06-12 10:33 | Batch | 2014-06-12 13:30 | Batch | 2014-06-12 13:30 | Batch | 2014-06-12 13:30 |

| Sample Name | 14.4441.dox-1 | Type | Sample | Vial | 62 | | | | |
|----------------|---------------|--------------|------------------|---------|---------------------|----------------------|---------------------|---------|------|
| Vol. [µl] | 5 | | Concent | Diox | | | | | |
| Compound | RT [min] | Conc [ng/ml] | TEQ Conc [ng/ml] | LOO | Upper Bound [ng/ml] | Medium Bound [ng/ml] | Lower Bound [ng/ml] | WIKATEF | 2005 |
| 2378-TCDD | 20.94 | -LOQ | -LOQ | 0.0200 | 0.02000 | 0.01000 | 0.00000 | 1 | |
| 12378-PCDD | 24.27 | -LOQ | -LOQ | 0.0300 | 0.03000 | 0.01500 | 0.00000 | 1 | |
| 123478-HCDD | 28.02 | 0.0376 | 0.00376 | 0.0200 | 0.00376 | 0.00376 | 0.00376 | 0.1 | |
| 123678-HCDD | 28.15 | -LOQ | -LOQ | 0.1400 | 0.01400 | 0.00700 | 0.00000 | 0.1 | |
| 123789-HCDD | 28.51 | -LOQ | -LOQ | 0.0600 | 0.00600 | 0.00300 | 0.00000 | 0.1 | |
| 1234678HpCDD | 33.08 | -LOQ | -LOQ | 0.4400 | 0.00440 | 0.00220 | 0.00000 | 0.01 | |
| OCDD | 39.46 | 4.9892 | 0.00150 | 1.9900 | 0.00150 | 0.00150 | 0.00150 | 0.0003 | |
| 2378-TCDF | 19.90 | -LOQ | -LOQ | 0.1000 | 0.01000 | 0.00500 | 0.00000 | 0.1 | |
| 12378-PCDF | 23.26 | -LOQ | -LOQ | 0.5600 | 0.01680 | 0.00840 | 0.00000 | 0.03 | |
| 23478-PCDF | 24.08 | -LOQ | -LOQ | 0.0800 | 0.02400 | 0.01200 | 0.00000 | 0.3 | |
| 123478-HCDF | 27.03 | -LOQ | -LOQ | 0.0700 | 0.00700 | 0.00350 | 0.00000 | 0.1 | |
| 123678-HCDF | 27.19 | -LOQ | -LOQ | 0.0600 | 0.00600 | 0.00300 | 0.00000 | 0.1 | |
| 123789-HCDF | 29.08 | -LOQ | -LOQ | 0.1800 | 0.01800 | 0.00900 | 0.00000 | 0.1 | |
| 234678-HCDF | 27.84 | -LOQ | -LOQ | 0.0800 | 0.00800 | 0.00400 | 0.00000 | 0.1 | |
| 1234678HpCDF | 31.29 | -LOQ | -LOQ | 0.5500 | 0.00550 | 0.00275 | 0.00000 | 0.01 | |
| 1234789HpCDF | 36.01 | -LOQ | -LOQ | 0.9200 | 0.00920 | 0.00460 | 0.00000 | 0.01 | |
| OCDF | 39.89 | -LOQ | -LOQ | 0.6000 | 0.00018 | 0.00009 | 0.00000 | 0.0003 | |
| PCB 77 | 18.01 | -LOQ | -LOQ | 49.6000 | 0.00497 | 0.00248 | 0.00000 | 0.0001 | |
| PCB 81 | 18.01 | 6.7691 | 0.00203 | 3.6700 | 0.00203 | 0.00203 | 0.00203 | 0.0003 | |
| PCB 126 | 21.04 | -LOQ | -LOQ | 1.3700 | 0.13700 | 0.06850 | 0.00000 | 0.1 | |
| PCB 169 | 24.22 | -LOQ | -LOQ | 0.0700 | 0.00210 | 0.00105 | 0.00000 | 0.03 | |
| Total TEQ | | | | | 0.32144 | 0.16436 | 0.00729 | | |
| Sum TEQ PCDD/F | | | | | 0.17534 | 0.09030 | 0.00526 | | |
| Sum TEQ PCBs | | | | | 0.14610 | 0.07406 | 0.00283 | | |

| Sum TEQ PCDD/F | 0.17534 |
|----------------|---------|
| Sum TEQ PCBs | 0.14610 |

Take Home Message

- ✓ PTV-GC/MS/MS accepted as a confirmatory tool under EU Regs
- ✓ Full validation on challenging matrix
- ✓ MS/MS, but still dioxin analyses...
- ✓ MS/MS & sectors to be properly perceived