

Investigation of Ion Mobility coupled with mass spectrometry (IMMS) for the screening of pesticide residues in food

June 2011, Mass Spectrometry in Food and Feed



S. Goscinny, R. Touilloux, L. Joly, E. De Pauw,
V. Hanot and G. Eppe




Our Path today:

- Purpose
- Power
- Pesticides
- Prospect


Our Path today:

- Purpose
- Power
- Pesticides
- Prospect

Multiresidue Methods are required tools for the determination of a great number of various compounds in one analysis



Has to be viable for the lab



Analytical trend for MRM

Sample
treatment



Separation



Measurement



Data
processing

Sample
treatment



Separation



Measurement



Data
processing

Generic
No purification



Dirty samples



Sample
treatment



Separation



Measurement



Data
processing

Generic
No purification



Fast
UPLC



Dirty samples



Co-elutions



Sample
treatment



Separation



Measurement



Data
processing

Generic
No purification



Fast
UPLC



Mass
spectrometry
revolution



Friendly
software

Dirty samples



Co-elutions




Optimization



Under control
No false + nor -

But...


Optimize the acquisition






But...

Optimize the acquisition



But...


Optimize the acquisition



A person is seen from behind, walking up a steep, rocky mountain trail. They are carrying two large, full white sacks on their back, secured with a dark strap. The person is wearing blue jeans and dark shoes. The trail is narrow and made of light-colored dirt and rocks. The surrounding environment is rugged, with steep, rocky slopes and sparse vegetation, including some green shrubs and bare trees. In the background, a steep mountain peak is visible under a clear sky. The text "Labor intensive" is overlaid in red on the upper part of the sacks.

**Labor
intensive**

**What's
next?**



Screening

Screening method


Method to **detect** analyte (s)
in samples
in an easy/rapid way

- Qualitative
- No full identification
- No value with uncertainty but result with a level of confidence


How can we help the
Detection ?


Our Path today:


- Purpose
- Power
- Pesticides
- Prospect




Synapt G2










Ion Mobility



Small, compact




Large extended



- A 80 years old concept (C.F. Powell, 1932)
- The idea: ions «race»; the most mobile reach the detector first
- Separation is driven by electric fields not under vacuum

Traditional IMS

Travelling Wave Ion Guides



Introduced by Waters with the range of Premier MS:

Eliminates crosstalk problems ➡ faster ion transit

T-Wave and IM separation

High electric field applied
SEQUENTIALLY through the IM
cell


Ions are moved through the
IM cell in **PULSES** as **WAVES**



faster IMS duty cycles




MS sensitivity is not
compromized



Low Mobility Ions

High Mobility Ions



Our Path today:

- Purpose
- Power
- Pesticides
- Prospect

Pesticides


Broad term for products that are used to control pests


Rich in diversity;


Chemical structure, solubility, volatility, potential
for degradation...





Perfect


Challenge









Compound	m/z	Structure
Methamidophos	141.0013	
Dichlorvos	219.9459	
Mepanipyrim	223.1109	
Spinosad	731.4608	





IM cell



Nature of the gas

IMS T-Wave velocity
(m/sec)




IMS T-Wave Height
(V)


Gas Pressure (mbar)

Biais (V)

Helium Cell Pressure
(mbar)

Nature of the Gas





Nature of the gas : N_2

IMS T-Wave velocity
(m/sec)


IMS T-Wave Height
(V)

Gas Pressure (mbar)

Biais (V)

Helium Cell Pressure

Last voltage before the IM cell



Bias (V)

Too low




bad transmission

Too high



Ion fragmentation


Buffer gas: softly reduces ion velocity



Maximises transmission of ions on entry into the IM cell

Helium Cell Pressure

IMS T-Wave velocity



Fast




Ions roll over the wave



Better the separation

IMS T-Wave Height



High field



High pulses




Better separation
selectivity

Gas Pressure

Higher



Better the resolving
power



Plackett-Burman design

So called « Screening Designs »

Finds influencing factors with a limited
number of experiments

Plackett-Burman design

5 parameters

IMS T-Wave velocity

IMS T-Wave Height

Gas Pressure

Helium Cell Pressure

Biais (V)

3 responses


Intensity

Resolution

Relative drift time

The construction
of the design is
done with 15 runs

Most influencing factors



Nature of the gas

IMS T-Wave velocity
(m/sec)


IMS T-Wave Height
(V)

Gas Pressure (mbar)

Biais (V)

Helium Cell Pressure

Most influencing factors



Gas Pressure (mbar)


Biais (V)


Optimization of 2 parameters!

Central Composite Design


The 3 less influencing parameters are set to the values of maximum separation


Then the CCD will be performed with the 2 most influencing parameters





Ion Mobility of 5 classes of pesticide






Nature's best

From 510 LC amenable pesticides, we only have 5 with a mass > 600


- 4 belong to the avermectins family (masses > 800 Da)
- Azadirachtin is isolated from Neem tree >700 Da)



Natural products/synthetic analogues



Matrix effect on the drift time




Interesting case study

Unusual Suspect


Quinalphos

299.0619 (+1)




Unusual Suspect

Quinalphos 299.0619 (+1)




Phenthoate 321.0384 (+1)




Unusual Suspect

Quinalphos 299.0619 (+1)




Phenthoate 321.0384 (+1)




Confirmation by MS/MS

Investigation

Quinalphos 299.0619 (+1)



Phenthoate 321.0384 (+1)



NO Phenthoate !




Investigation

Quinalphos

299.0619 (+1)


8.44
(min)



Phenthoate

321.0384 (+1)

8.48
(min)




Co-elution


<p>11 Na Sodium 22.989770 [Ne]3s 5.1391</p>	<p>12 Mg Magnesium 24.3050 [Ne]3s² 7.6462</p>	<p>20 Ca Calcium 40.078</p>
<p>19 K Potassium 39.0983</p>	<p>21 Sc Scandium 44.9559</p>	

Unusual Suspect

Quinalphos 299.0619 (+1)




Phenthoate 321.0384 (+1)




Quinalphos
+ Na 321.0439 (+1)

Unusual Suspect

Quinalphos 299.0619 (+1)




Phenthoate 321.0384 (+1)




Quinalphos
+ Na 321.0439 (+1)

Can IMS help?

Unusual Suspect



Unusual Suspect



Our Path today:

- Purpose
- Power
- Pesticides
- Prospect

Future

- Many optimization options, many possibilities
- IM adds value to the data obtained, ...new IP?
- Need to test the UPLC-IMMS application
- Investigate the impact of the chemical functions of the different pesticides on the drift behavior

... and let our curiosity

and imagination run free!!

