

Reading the decomposition chemistry of human remains with a new pair of glasses

Stefanuto PH^a, Stadler S^b, Perrault K^c, Pesesse R^a, Forbes S^c, Focant JF^a

^aOrganic and Biological Analytical Chemistry, Chemistry Department, University of Liège, Allée de la Chimie 3, B-6c Sart-Tilman, B-4000 Liège, Belgium.

^bFaculty of Science, University of Ontario Institute of Technology, 2000 Simcoe St N., Oshawa, Ontario, Canada L1H 7K4.

^cCentre for Forensic Science, University of Technology, Sydney, PO Box 123, Broadway, NSW, 2007, Australia.

For the last few decades, investigations have been conducted to better learn the process of body decomposition by mean of the measurement of the VOCs released during decay. However, the chemical fingerprint of the decomposition odor is still far from being elucidated. We recently developed an analytical strategy that relies on the use of thermal desorption (TD) coupled to comprehensive two-dimensional GC coupled to time-of-flight MS (GC×GC-TOFMS) for the profiling of cadaveric VOCs from grave soils and decaying bodies.^{1,2} The additional peak capacity of GC×GC, the spectral deconvolution algorithms applied to unskewed mass spectral data, and the use of robust specific data mining strategies allowed the generation of characteristic VOC profiles across the various stages of soft-tissue decomposition.³

For practical and ethical reasons, human analogs (*Sus domesticus* L. carcasses) were used in these early investigations. For the present study, we investigated the VOC profile of early stage decomposition of human bodies. We analyzed samples collected during different trials organized during different seasons in a body farm located in Texas. Samples included environmental controls, pig carcasses, and human bodies (protected or not from scavenger insects). This approach is a step forward in producing specific complex odor profiles that could be used in various types of forensic investigations.

Human remains detection (HRD) canines utilize this odor signature to locate human remains during police investigations and recovery missions in the event of a mass disaster. Currently, it is unknown what compounds or combinations of compounds are recognized by the HRD canines.⁴ The better understanding of the VOC profile released by dead or injured bodies will possibly help to designing more realistic canine artificial scent training solutions, but also allow enhancing the specificity of electronic noses which, so far, still suffers from poor efficiency at least partly due to the lack of information about the decomposition chemistry.

¹Dekeirsschieter, J. et al. PLoS ONE (2012) 7, e39005.

²Brasseur C. et al. J Chromatogr A (2012) 1255, 163–170.

³Stadler, S. et al. Anal Chem (2013) 85, 998–1005.

⁴Stadler, S. et al. J Chromatogr A (2012) 1255, 202–206.