

A fast gas chromatographic method for the study of semiochemical slow release formulations

S. Heuskin¹, S. Lorge², P. Leroy³, F. Verheggen³, J.-P. Wathelet⁴, G. Lognay¹

¹ Department of Analytical Chemistry, University of Liege, Gembloux Agro-Bio Tech, Passage des Déportés 2, B-5030 Gembloux, Belgium

² Department of Interfaces Chemistry, UCL, Croix du Sud 1, B-1348 Louvain-la-Neuve, Belgium

³ Department of Functional and Evolutive Entomology, University of Liege, Gembloux Agro-Bio Tech, Passage des Déportés 2, B-5030 Gembloux, Belgium

⁴ Department of General and Organic Chemistry, University of Liege, Gembloux Agro-Bio Tech, Passage des Déportés 2, B-5030 Gembloux, Belgium

Introduction


E- β -Farnesene, the alarm pheromone of many aphid species and β -caryophyllene, recently identified as one of the possible component of the aggregation pheromones of the Asian ladybeetles *Harmonia axyridis* Pallas, are considered as two sesquiterpenes attractive for aphids' predators and parasitoids.

In the present research, a fast GC method was optimised for the quantitative characterisation, in less than five minutes, of the above mentioned semiochemicals formulated in slow release devices. The formulations were evaluated in terms of volatiles release capacity and protection efficiency of sesquiterpenes. The sesquiterpenes used in the formulations were obtained from natural sources. They were purified by flash chromatography fractionation (under pressure) of essential oils.

The fast chromatographic method was completely validated by means of the accuracy profile concept.

E- β -farnesene

β -caryophyllene



Experimental

Flash Chromatography fractionation of essential oils




Matricaria chamomilla

Compounds	
Sum of monoterpenes	1.3 %
E-β-farnesene	83.8 %
Germacrene D	1.5 %
Bicyclogermacrene	1.5 %
(E,E)- α -farnesene	11.9 %

Nepeta cataria

Compounds	
Sum of monoterpenes	1.2 %
β-caryophyllene	97.7 %
α -humulene	1.1 %

Fast GC analytical method



Calculated selectivity factor (α) between longifolene and β -caryophyllene : **1.016 > 1**

Chromatographic conditions

Ultra Fast Module : Ph5; 0.1 μ m film thickness, 5m x 0.1mm I.D.

Carrier gas : He; 0.5 ml/min

Split ratio : 1:100

Oven :

Initial T° : 40°C; 0.10 min

Ramp 1 : 30°C/min \rightarrow 95°C

Ramp 2 : 35°C/min \rightarrow 155°C


Ramp 3 : 200°C/min \rightarrow 280°C; 0.5 min

Oven run time : **4.78 min**

Validation by accuracy profile

Calibration standards


Linear regression models



E- β -farnesene


Validation standards

Back-calculation of the results \rightarrow linearity curves






Accuracy Profile

ACCURACY = TRUENESS + PRECISION





β -caryophyllene



Applications



Volatile release capacity of formulations



Laboratory controlled conditions :

- Temperature : 20°C
- Sampling air flow : 0.5 L/min
- Relative humidity : 65 %

Sesquiterpenes protection efficiency of formulations



- - - Alginate gel beads without alpha-tocopherol
- - - Alginate gel beads with alpha-tocopherol
- - - Sunflower oil formulation
- - - Purified compounds non formulated

Conclusions

- Fast achievement of high purity semiochemicals by flash chromatography of essential oils.
- Fast GC method for characterization and quantification of sesquiterpenes in less than 5 minutes for various applications. Moreover the good analytical performances of the fast column were demonstrated (see Heuskin et al., 2009, J. of Chromatography A, 1216, pp 2768-2775).
- Method completely validated for the two sesquiterpenes by the accuracy profile concept considering a risk of 5 % and a linear regression model.
- Low limits of quantification